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
Vol. 9

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- 1) 三ッ本祐樹: 超原子価ヨウ素反応剤とその利用. 化学と工業, 65; 330-331, 2012.
- 2) Green WA: *Industrial Photoinitiators*. CRC Press: Boca Raton, 2010.
- 3) Crivello JV, Lee JL: Alkoxy-substituted diaryliodonium salt cationic photoinitiators. *J. Polym. Sci., Part A: Polym. Chem.*, 27; 3951-3968, 1989.
- 4) Eick JD, Kostoryz EL, Rozzi SM, Jacobs DW, Oxman JD, Chappelow CC, Glaros AG, Yourtee DM: In vitro biocompatibility of oxirane/polyol dental composites with promising physical properties. *Dent Mater.*, 18; 413-421, 2002.
- 5) Crivello JV: The Discovery and development of onium salt cationic photoinitiators. *J. Polym. Sci., A: Polym. Chem.*, 37; 4241-4254, 1999.
- 6) Chen YC, Ferracane JL, Prahl SA: Quantum yield of conversion of the photoinitiator camphorquinone. *Dent. Mater.*, 23; 655-664, 2007.
- 7) Cook WD, Chen F: Enhanced photopolymerization of dimethacrylates with ketones, amines, and iodonium salts: The CQ system. *J. Polym. Sci., Part A: Polym. Chem.*, 49; 5030-5041, 2011.
- 8) Cook WD: Photopolymerization kinetics of dimethacrylates using the camphorquinone/amine initiator system. *Polymer*, 33; 600-609, 1992.
- 9) Jakubiak J, Allonas X, Fouassier JP, Sionkowska A, Andrzejewska E, Linden LÅ, Rabek JF: Camphorquinone-amines photoinitiating systems for the initiation of free radical polymerization *Polymer*, 44; 5219-5226, 2003.
- 10) Weinmann W, Thalacker C, Guggenberger R: Siloranes in dental composites. *Dent. Mater.*, 21; 68-74, 2005.
- 11) Shin DH, Rawls HR: Degree of conversion and color stability of the light curing resin with new photoinitiator systems. *Dent. Mater.*, 25; 1030-1038, 2009.
- 12) de Andrade KMG, Palialola AR, Lancellotti AC, Aguiar FHB, Watts DC, Gonçalves LS, Lima AF, Marchi GM: Effect of diphenyliodonium hexafluorophosphate on resin cements containing different concentrations of ethyl 4-(dimethylamino)benzoate and 2-(dimethylamino)ethyl methacrylate as co-initiators. *Dent. Mater.*, 32; 749-755, 2016.
- 13) Gonçalvesa LS, Moraes RR, Ogliari FA, Boaro L, Braga RR, Consani S: Improved polymerization efficiency of methacrylate-based cements containing an iodonium salt. *Dent. Mater.*, 29; 1251-1255, 2013.
- 14) Ogliari FA, Ely C, Petzhold CL, Demarco FF, Piva E: Onium salt improves the polymerization kinetics in an experimental dental adhesive resin. *J. Dent.*, 35; 583-587, 2007.
- 15) Dressano D, Palialol AR, Xavier TA, Braga RR, Oxman JD, Watts DC, Marchi GM, Lima AF: Effect of diphenyliodonium hexafluorophosphate on the physical and chemical properties of ethanolic solvated resins containing camphorquinone and 1-phenyl-1,2-propanedione sensitizers as initiators. *Dent. Mater.*, 32; 756-764, 2016.
- 16) Ogliari FA, Ely C, Lima GS, Conde MCM, Petzhold CL, Demarco FF, Piva E: Onium salt reduces the inhibitory polymerization effect from an organic solvent in a model dental adhesive resin. *J. Biomed. Mater. Res. Part B: Appl. Biomater.*, 86B; 113-118, 2008.
- 17) Leal FB, Lim GS, Collares FM, Samuel SM, Petzhold CL, Piva E, Ogliari FA: Iodonium salt improves the dentin bonding performance in an experimental dental adhesive resin. *Inter. J. Adhes. Adhes.*, 38; 1-4, 2012.
- 18) Kim M, Suh B, Shin D, Kim K-M: Comparison of the physical and mechanical properties of resin matrix with two photoinitiator systems in dental adhesives. *Polymers*, 8; 250, 2016.
- 19) Abedin F, Ye Q, Song L, Ge X, Camarda K, Spencer P: Effect of partition of photo-initiator components and addition of iodonium salt on the photopolymerization of phase-separated dental adhesive. *JOM.*, 68; 1090-1099, 2016.
- 20) Guo X, Peng Z, Spencer P, Wang Y: Effect of initiator on photopolymerization of acidic, aqueous dental model adhesives. *J. Biomed. Mater. Res.*, 90A; 1120-1127, 2009.
- 21) Wang Y, Spencer P, Yao X, Ye Q: Effect of coinitiator and water on the photoreactivity and photopolymerization of HEMA/camphorquinone-based reactant mixtures. *J. Biomed. Mater. Res.*, 78A; 721-728, 2006.
- 22) Oguri M, Yoshida Y, Yoshihara K, Miyauchi T, Nakamura Y, Shimoda S, Hanabusa M, Momoi Y, Van Meerbeek B: Effects of functional monomers and photo-initiators on the degree of conversion of a dental adhesive. *Acta Biomater.*, 8; 1928-1934, 2012.
- 23) Hanabusa M, Yoshihara K, Yoshida Y, Okihara T, Yamamoto T, Momoi Y, Van Meerbeek B: Interference of functional monomers with polymerization efficiency of adhesives. *Eur. J. Oral Sci.*, 124; 204-209, 2016.
- 24) Salz U, Mücke A, Zimmermann J, Tay FR, Pashley DH: pKa value and buffering capacity of acidic monomers commonly used in self-etching primers. *J. Adhes. Dent.*, 8; 143-150, 2006.
- 25) Tay FR, Pashley D, Suh BI, Hiraishi N, Yiu CKY: Water treeing in simplified dentin adhesives-Déjà Vu? *Oper. Dent.*, 30; 561-579, 2005.
- 26) Loguercio AD, Stanislawczuk R, Mittelstadt FG, Meier MM, Reis A: Effects of diphenyliodonium salt addition on the adhesive and mechanical properties of an experimental adhesive. *J. Dent.*, 41; 653-658, 2013.
- 27) Besse V, Pluart LL, Cook WD, Pham T-N, Madec P-J: Polymerization kinetics of phosphonic acids and esters using an iodonium initiator. *J. Polym. Sci., Part A: Polym. Chem.*, 51; 5046-5055, 2013.
- 28) Feitosa VP, Sauro S, Ogliari FA, Stansbury JW, Carpenter GH, Watson TF, Sinhoreti MA, Correr AB: The role of spacer carbon chain in acidic functional monomers on the physicochemical properties of self-etch dental adhesives. *J. Dent.*, 42; 565-574, 2014.
- 29) Meereis CTW, Leal FB, Ogliari FA: Stability of initiation systems in acidic photopolymerizable dental material. *Dent. Mater.*, 32; 889-898, 2016.
- 30) Guo X, Wang Y, Spencer P, Ye Q, Yao X: Effects of water content and initiator composition on photopolymerization of a model BisGMA/HEMA resin. *Dent. Mater.*, 24; 824-831, 2008.
- 31) Ullrich G, Burtcher P, Salz U, Moszner N, Liska R: Phenylglycine derivatives as coinitiators for the radical photopolymerization of acidic aqueous formulations. *J. Polym. Sci., Part A: Polym. Chem.*, 44; 115-125, 2006.
- 32) Morlet-Savary F, Klee JE, Pfefferkorn F, Fouassier JP, Lalevée J: The camphorquinone/amine and camphorquinone/amine/phosphine oxide derivative photoinitiating systems: Overview, mechanistic approach, and role of the excitation light source. *Macromol. Chem. Phys.*, 216; 2161-2170, 2015.

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1. Yamada B, Zetterlund PB, “General chemistry of radical polymerization”. In *Handbook of Radical Polymerization*, Matyjasewski K, Davis TP, eds, Wiley-Interscience, New York, pp.117-186, 2002.
2. Yamada B, Zetterlund PB, Sato E, Utility of propenyl groups in free radical polymerization: Effects of steric hindrance on formation and reaction behavior as versatile intermediates. *Prog. Polym. Sci.*, 31:835-877, 2006.
3. Kubota B, Kajiwara A, Zetterlund PB, Kamachi M, Treurnicht J, Tonge MP, Gilbert RG, Yamada B: Determination of the propagating radical concentration. *Macromol. Chem. Phys.*, 208: 2403-2411, 2007.
4. McHale R, Carroll WM, Aldabbagh F, Yamada B, A Study of efficient synthesis and copolymerization of polyacrylic acid and polyacrylic ester macromonomers: Manipulation by steric factor. *Macromol. Chem. Phys.*, 206: 2054-2066, 2005.

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